


# FluxMed

## Respiratory Mechanics Monitor

Measurement in all the phases of ventilation




### MAKES POSSIBLE

- Faster weaning
- Shorter ventilation time
- Shorter hospitalization time
- Greater successful weaning rate
- Avoid complications associated with mechanical ventilation
- Optimize the patient-ventilator synchrony
- Cost savings


### ADVANCED MEASUREMENT

- Transpulmonary Pressure
- Work of Breathing (WOB)
- Campbell Diagram
- Pressure Time Product (PTP)
- Stress Index (SI)
- Transpulmonary Stress Index
- Volumetric capnography
- Alveolar CO<sub>2</sub> (PACO<sub>2</sub>)
- Breath by breath Vd/Vt
- CO<sub>2</sub> production
- Percentage of assistance
- Electrical Impedance Tomography



## SPECIALLY DESIGNED TO MEASURE


Controlled mandatory ventilation  
Assist-control ventilation  
Weaning process  
T-Piece Test  
Non Invasive Ventilation  
Spontaneous Breathing  
Rehabilitation and Follow-Up



## WORK OF BREATHING

Pressure Time Product (PTP) makes possible to measure the patient respiratory effort, divided in the elastic and resistive components.

The Campbell diagram shows the patient work of breathing (WoB), separating the different elastic and resistive areas. These tools are especially useful for complex weaning and for COPD patients.



## STRESS INDEX

Stress Index indicates if a patient is being ventilated in the safe zone, avoiding overdistension and alveolar collapse. Stress Index is especially useful for ARDS patients.

## VOLUMETRIC CAPNOGRAPHY

Makes possible to assess the efficiency of the patient's ventilation. By knowing the dead space to tidal volume ratio ( $V_d/V_t$ ), the  $CO_2$  alveolar pressure and the  $CO_2$  production ( $VCO_2$ ), breath by breath information about patient's lung perfusion can be obtained.



## ELECTRICAL IMPEDANCE TOMOGRAPHY

Adding EIT to the respiratory monitoring provides a deeper understanding of the physiology of patients requiring mechanical ventilation. This non-invasive technology allows a regional and continuous analysis of the lungs. Fluxmed EIT uses SenTec AG EIT Branch technology.

## SIGNALS AND PARAMETERS RECORDING AND REVIEWING

The FluxView software allows to record and view respiratory mechanics signals and parameters in real time. The FluxReview software makes possible to analyze previously recorded signals.