


# Disposa-View™



The **ONLY COMPLETE** Disposable  
Head Positioning Device!



Designed for the  
Prone-Positioned Patient

vital **signs**<sup>®</sup>

# Protect Your Prone-Positioned Patients

A latex-free single-use head positioning device for prone procedures.



The soft, comfortable foam reduces the pressure around the patient's eyes and face, while supporting the head and neck.

The mirror provides the clinician with a clear view of the face, eyes, and endotracheal tube for easy monitoring.

The slotted design allows for the endotracheal tube to be positioned on either side of the patient's face.



Product Specifications and Ordering Information	
Product Name	Disposa-View™
Product Code	8000HDP
Description	Prone Position Head Cushion with Mirror
Quantity	6/Case

20 Campus Road, Totowa, NJ 07512, USA • TEL:800-932-0760 • FAX:+1-973-790-3307 • www.vital-signs.com  
Authorized European Representative: Vital Signs, Ltd., 13-14 Eidon Way, Lineside Industrial Estate, Littlehampton, West Sussex, UK

