



CARETEK

Pacific Couches

PE320



PE120



PE220






Features:

- Vertical Lift (no horizontal displacement)
- Tente 5" Central Locking Castors as Standard
- Mobile Hoist Compatible
- Optional - Fold Down Side Supports
- Optional - Paper Roll Holder

Length	Width	Height Adjustment Range (Linak Electric Adjustment)
1950mm	700mm	460mm to 950mm
Back Rest Adjustment (Linak Electric Adjustment)	Rest Adjustment (Linak Electric Adjustment)	Safe Working Load
0° to 80°	-30° to 0°	250kg



Parameter comparison



PE320



PE220



PE120

	PE320	PE220	PE120
Backrest Angle Adjustment (Electric)	●	●	
Footrest Angle Adjustment (Electric)	●		
Vertical Lift (No Horizontal Displacement) (Electric)	●	●	●
ABS Plastic Undercarriage Cover	●	●	●

Standard ●