


CBLGROUP


Caretek medical & care
over bed tables

LEAN SIX-SIGMA





Customers



ARJOHUNTLEIGH
GETINGE GROUP



NHG
NURSING HYGIENE GROUP
CARING ABOUT CARE



CBLGROUP

- Flat packed
- Smooth, mechanical spring assisted height adjustment
- Positive lock at any selected height
- Fixed one piece and two piece part tilting table options
- Tilting top, adjustable to a variety of angles thanks to the light weight, easy to use Hettich ratchet system
- Smooth non absorbent tabletop surfaces for easy cleaning
- Safe working load 25 kg evenly distributed
- 12 months warranty / CE Marked





MT110 Fixed Top Table



- Mild steel construction powder coated off white
- Vacuum laminated MDF timber top light blue or white
- Top size 900 x 425 (mm)
- Locks in compression no release required for raising
- Modern ABS plastic base cover
- 55mm Diameter castors
- Adjustment range 745 to 1055 (+/- 5mm)
- Safe working load 25kg evenly distributed

Loading Capacity :

40' container - 912 pieces

40' HC container - 990 pieces





MT210 Split Top Table



- Mild steel construction powder coated off white
- Vacuum laminated MDF timber top light blue or white
- Top size 700 x 425 and 195 x 400 (mm)
- Locks in compression no release required for raising
- Modern ABS plastic base cover
- 55mm Diameter castors
- Tilt Section angles 24,31,37,42,47 degrees
- Adjustment range 745 to 1055 (+/-5mm)
- Safe working load 25kg evenly distributed
- Hettich bracket

Loading Capacity:

40' container - 912 pieces

40' HC container - 990 pieces





MT210 - BT Split Top Table

- Mild steel construction powder coated off white
- Laminated beech MDF timber top with melamine veneer
- Top size 700 x 425 and 195 x 400 (mm)
- Locks in compression no release required for raising
- Modern ABS plastic base cover
- 55mm Diameter castors
- Tilt Section angles 24,31,37,42,47 degrees
- Adjustment range 745 to 1055 (+/-5mm)
- Safe working load 25kg evenly distributed
- Hettich bracket

PU finish / Standard matt

Loading Capacity:

40' container - 912 pieces

40' HC container - 990 pieces





KT110 – N2 Fixed Top Table



Loading Capacity:

40' container - 912 pieces

40' HC container - 990 pieces

- Mild steel construction powder coated off white
- Laminated beech MDF timber top with melamine veneer. 2 strips
- Top size 885 x 395 (mm)
- Locks in compression no release required for raising
- Modern ABS plastic base cover
- 55mm Diameter castors
- Adjustment range 738 to 1048 (+/- 5mm)
- Safe working load 25kg evenly distributed

PU finish / Standard matt





Loading Capacity:

40' container - 912 pieces

40' HC container - 990 pieces

HT110-N2 Fixed Top Table

- Mild steel construction powder coated off white
- Laminated beech MDF timber top includes 2 melamine veneer strips
- Top size 885 x 395 (mm)
- Locks in compression no release required for raising
- 55mm Diameter castors
- Adjustment range 760 to 1070 (+/- 5mm)
- Safe working load 25kg evenly distributed

PU gloss finish / standard matt



HT110-P Fixed Top Table

- Mild steel construction powder coated off white
- Blow moulding table top
- Top size 885 x 395 (mm)
- Locks in compression no release required for raising
- 55mm Diameter castors
- Adjustment range 760 to 1070 (+/- 5mm)
- Safe working load 25kg evenly distributed

Loading Capacity:

40' container - 912 pieces

40' HC container - 990 pieces



Loading Capacity:

40' container - 912 pieces

40' HC container - 990 pieces

HT110-S2 Fixed Top Table (short version)

- Mild steel construction powder coated off white
- Laminated beech MDF timber top includes 2 melamine veneer strips
- Top size 885 x 395 (mm)
- Locks in compression no release required for raising
- 55mm Diameter castors
- Adjustment range 630 to 825 (+/- 5mm)
- Safe working load 25kg evenly distributed

PU gloss finish / Standard matt

BT110 Fixed Top Table



Loading Capacity:

40' container - 912 pieces

40' HC container - 990 pieces

- Mild steel construction powder coated off white/grey
- Laminated beech MDF timber top
- Top size 885 x 395 (mm)
- compression spring in inside, locked by hand wheel. Release the hand wheel if require raising
- 55mm Diameter castors
- Adjustment range 630 to 825 (+/- 5mm)
- Safe working load 25kg evenly distributed

PU gloss finish / Standard matt



2021

CERTIFICATE OF REGISTRATION

This certifies that:

Caretek(China) Medical PLC
Xiguakeng, Guanqiao Village, Shilou Town
Panyu District
Guangzhou City Guangdong, CN 511447

is registered with the U.S. Food and Drug Administration for FY 2021 pursuant to Title 21, 807 et seq. of the United States Code of Federal Regulations:


Establishment Registration: **3010509576**
Device Classification Name: **CHAIR AND TABLE, MEDICAL**
Product Code: **KMN**
Regulation Number: **880.6140**
Official Correspondent and U.S. Agent: **Registrar Corp**
144 Research Drive, Hampton, Virginia, 23666, USA
Telephone: +1-757-224-0177 • Fax: +1-757-224-0179

Registrar Corp will confirm that such registration remains effective upon request and presentation of this certificate until the end of the year stated above, unless said registration is terminated after issuance of this certificate. Registrar Corp makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. Registrar Corp assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding."

The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration. Registrar Corp is not affiliated with the U.S. Food and Drug Administration.


Registrar Corp
144 Research Drive, Hampton, Virginia, 23666, USA
Telephone: +1-757-224-0177 • Fax: +1-757-224-0179
info@registrarcorp.com • www.registrarcorp.com


David Lemarz
Executive Director
Registrar Corp
Dated: *October 28, 2020*

©2005-2020 Registrar Corp

Certificate – 證明書 – Сертификат – Certificat – 증명서 – certificado

CERTIFICATE OF CONFORMITY



Certificate No. CTK20160927002GPSD

Certificate's Holder : **Caretek (China) Medical PLC**
Address : Xiguakeng, Guanqiao Village, Shilou Town,
Panyu District, Guangzhou, Guangdong Province,
P.R.China

Product Name : **Overbed Table**
Model(s) : MT210, MT211, MT110, MT111, MT112, KT110, KT111,
HT110, HT111, ST110, ST111, CT110, CT111, TT110,
TT111, HT110-AC-01, HT110-MT, HT110-S2, HT110P,
KT210

Applicable Standard(s) : EN 581-1:2017
EN 581-3:2017
BS 2483:1997

The submitted products have been tested by us with the above standards and found in compliance with the following European Directives:

General Product Safety Directive 2001/95/EC

Note: This certificate is only valid for the equipment and configuration described, and in conjunction with the test data detailed above. The CE mark as show joined can be used, under the responsibility of the manufacturer or the importer, after completion of the CE Declaration of Conformity and in accordance with the above directives.

Issuing Date: 06.11.2019
Expiry Date: 26.09.2021



Thomas Diong
Thomas Diong
Chief Certification Manager

GUANGZHOU CERTITEK TESTING SERVICES CO., LTD.

Room703/7F, Development Building, Tian An Hi-Tech Ecological Park

No.555 North Road Panyu Avenue, Panyu District, Guangzhou, 511400, China

+86 20 3921 1670 +86 20 3921 1640 info@certitek.cn www.certitek.co.uk v3.0