



Australia Caretek Medical

Our products

Home care
range--Table
and Cabinet



Product Code	Description
LP110	Bedside Locker



Overall size – W465 x D465 x H865 (mm)

- ABS Plastic construction
- Colors, off white & green or blue
- Side panels include drop down towel rail
- Pull out drawer and left hand cupboard door
- Open bottom shoe shelf



** Pull out drawer can be turned to create a small writing board



*****Two size*****

Product Code	Description	US\$
LPS100	Bedside Locker	\$ 135



LpS100 W565 x D495 x H755 (mm)

Mild steel construction

Colours, beige or off white

Top is high quality laminated MDF in beech or grey

Pull out drawer and left or right hand cupboard door

Product Code	Description	
L100	Bedside Locker	



Colours, light blue top, grey infill panels
 Top is made from quality MDF with vacuum form laminating
 Side panels are made from laminated MDF
 Aluminium corner extrusions
 Base is a welded mild steel, powder coated construction
 Drawer face laminated MDF.
 Internal drawer injection moulded



Overall size - W500 x D495 x H870 (mm)