



Compact Regulator with Yoke



Compact Regulator with Nut & Nipple



Click-Style Regulator with Yoke



Click-Style Regulator with Nut & Nipple



# Gas Regulators

## Click Style and Compact Gas Regulators

### Product Benefits

- Fully customizable to meet your facility's specific needs
- Customize Regulator style, connection and gas
- Choose from a variety of colors
- ISO or USA option
- Built in Flowmeter in click style option
- Add traditional Thorpe Style Flowmeter to Compact Regulator if required
- Meets all CGA, HCA, FDA, CE and ISO Requirements
- 100% tested and quality inspected
- 5 year warranty



Visit [www.amvex.com](http://www.amvex.com)  
for the Gas Regulator part  
configurator and much more!


## Product Features

### Click Style



- ① Amvex Click Style is available in 50 PSI Preset
- ② Available with a CGA Nut and Nipple or CGA Yoke Style
- ③ Available in Oxygen only
- ④ Aluminum body with brass core to ensure safety
- ⑤ Built in Flowmeter available in 0-4 LPM, 0-8 LPM, 0-15 LPM, and 0-25 LPM increments
- ⑥ Available in a variety of colors (blue is standard)

### Compact Style



- ① Amvex Compact Regulators are available in 50 PSI Preset
- ② Available with a CGA Nut and Nipple or CGA Yoke Style
- ③ Available in Oxygen, Air, N<sub>2</sub>O, and CO<sub>2</sub>
- ④ Available in USA and ISO colors
- ⑤ Available with standard Thorpe Style Flowmeter options
- ⑥ Machined brass



## Part Configurator

### Click Style Regulators

#### GRXX-YYY-(Z)

| Gas Regulator Type:<br>Flow Range (litres per minute) |    |
|---|----|
| 4 LPM:  | 4  |
| 8 LPM:  | 8  |
| 15 LPM:   | 15 |
| 25 LPM:   | 25 |

| Connection Type:  |     |
|-------------------|-----|
| CGA Yoke:         | 870 |
| CGA Nut & Nipple: | 540 |

| Optional Colors:        |   |
|-------------------------|---|
| Purple:                 | P |
| Green:                  | G |
| Gold:                   | Y |
| Red:                    | R |
| Black:                  | B |
| Silver:                 | S |
| <i>Blue is standard</i> |   |

### Compact Regulators

#### GR-XXX-(YYYY)

| Gas and Connection Type:       |     |
|--------------------------------|-----|
| Oxy Nut & Nipple:              | 540 |
| Oxy Yoke Style:                | 870 |
| Air Nut & Nipple:              | 346 |
| Air Yoke Style:                | 950 |
| N <sub>2</sub> O Nut & Nipple: | 326 |
| N <sub>2</sub> O Yoke Style:   | 910 |
| CO <sub>2</sub> Nut & Nipple:  | 320 |
| CO <sub>2</sub> Yoke Style:    | 940 |

| Optional Flowmeter: |       |
|---------------------|-------|
| 15 LPM Flowmeter:   | FM15  |
| 8 LPM Flowmeter:    | FM08  |
| 3.5 LPM Flowmeter:  | FM03  |
| 1 LPM Flowmeter:    | FM01  |
| 200 CCM Flowmeter:  | FM200 |