

# Springfield Crib

Made in USA

October 2013



Springfield Crib with optional IV Pump Bracket

## Features

- 3" Foam Mattress (CA Code TB 603)
- LAGX Gatch Extender
- E-Z Lift Sides
- 5" Casters (Set of 4, 2 with brakes)
- Clear Vinyl Panels (4 sides)
- Confetti Epoxy Finish
- Anti-Microbial Coating
- Steel Construction

## Options

- IV Pump Bracket
- Undercrib Shelf
- 4" Foam Mattress
- 5" Innerspring Mattress
- Vinyl Top Cover

## Articulation:

Head/Thigh/Foot/Section                      35/35/20 degrees

### HEMOCARE SPRINGFIELD CRIB

MODEL #	H1943	H1955	H1975
Crib Dimensions - Overall	30" x 60"	36" x 60"	36" x 72"
Mattress Size	28" x 58"	34" x 58"	34" x 70"
Overall Crib Height	75"	75"	75"
Floor to Top of Platform	29"	29"	29"
Height of Sides	45"	45"	45"

All Springfield cribs have a 2 year warranty.  
 HARD Manufacturing Co., Inc. is ISO 13485 and ISO 9001 certified

For additional information call 1-800-873-4273 or e-mail: [hardmfg@hardmfg.com](mailto:hardmfg@hardmfg.com)

