

Standard Crib

Made in USA

October 2013



Features

- 3" Foam Mattress (CA Code TB 603)
- E-Z Lift Sides
- 5" Casters (Set of 4, 2 with brakes)
- Confetti Epoxy Finish
- Anti-Microbial Coating
- Steel Construction

Options

- IV Pump Bracket
- Undercrib Shelf
- 4" Foam Mattress
- 5" Innerspring Mattress
- Pressure Reduction Mattress

Articulation:

Head/Thigh/Foot Section 35/35/20 degrees

<i>Homecare Standard Crib</i>				
MODEL #	H1942	H1952	H1973	H1982
Crib Dimensions - Overall	30" x 60"	36" x 60"	36" x 72"	36" x 83"
Mattress Size	28" x 58"	34" x 58"	34" x 70"	34" x 79"
Overall Crib Height	56"	56"	56"	56"
Floor to Top of Platform	31"	31"	31"	31"
Height of sides	26"	26"	26"	26"

All Standard cribs have a 2 year warranty.
 HARD Manufacturing Co., Inc. is ISO 13485 and ISO 9001 certified

For additional information call 1-800-873-4273 or e-mail: hardmfg@hardmfg.com

