



### Non-Rebreathing Mask

Delivering a high concentration of oxygen to spontaneously breathing patients, even with irregular breathing patterns and tidal volumes.

#### Reduced Risk of Patient Trauma

Elasticated earloop option on the adult mask eliminates the need for the patient's head to be moved to fit the mask, reducing risk of further injury.

#### Secured Oxygen Delivery

Star Lumen Tubing maintains oxygen flow, even if the tube is kinked.

#### Ease and Speed of Fitting

Mask is fitted quickly and easily with the swivel connector rotated to suit the patient's position.

#### Reduced Risk of Bag Occlusion

Castellations on the swivel connector ensures flow of oxygen from bag to patient is uninterrupted.

MASK	+ 2.1M TUBING		STRAP TYPE	QTY
032-10-071	032-10-070	Adult Non-Rebreathing Mask	Earloops	50
032-10-073	032-10-072	Adult Non-Rebreathing Mask	Head Strap	50
	032-10-072NC	Adult Non-Rebreathing Mask with Nose Clip	Head Strap	50
	032-10-0725M	Small Adult Non-Rebreathing Mask	Head Strap	50
	032-10-076	Adult Non-Rebreathing Mask with Safety Port	Earloops	50
	032-10-078	Adult Non-Rebreathing Mask with Safety Port	Head Strap	50
	032-10-078NC	Adult Non-Rebreathing Mask with Safety Port and Nose Clip	Head Strap	50
032-10-081	032-10-080	Paediatric Non-Rebreathing Mask	Head Strap	50
	032-10-084	Paediatric Non-Rebreathing Mask with Safety Port	Head Strap	50

### Tracheostomy Mask

Contoured for an effective seal and a 360° swivel connector for maximum patient comfort. The 22mm connector fits to corrugated tubing and can deliver humidified oxygen when used with a venturi and humidity adaptor (5506).



QTY

032-10-091	Adult Tracheostomy Mask	50
032-10-092	Paediatric Tracheostomy Mask	50

### Gilston T-Piece

For oxygen delivery to patients fitted with a tracheostomy tube. Sliding cover gives easy access for suction and a venturi barrel fitted to the tube to provide oxygen at a fixed concentration, or combined with a humidity adaptor for humidification.



QTY

5060	Adult Gilston T-Piece	5
5061	Paediatric Gilston T-Piece	5