



FCOT FIXED OXYGEN DELIVERY
The complete Post-op Solution.

FCOT is a complete solution for delivery of 40% oxygen to post-op adult and paediatric patients with a laryngeal mask, through recovery to the ward.

Visual Confirmation of Spontaneous Breathing

Condensation on the T-Piece and Indicator Bag inflation/deflation cycle clearly shows the patient's breathing pattern.

No Rebreathing or Oxygen Dilution

The high flow of oxygen enriched air inflates the Indicator Bag and flushes expired gases out through the vents, eliminating rebreathing. During inspiration the Indicator Bag acts as a reservoir, preventing dilution of the inhaled gas by ambient air.

Cost Effective

Complete kit with 40% Venturi, T-Piece with Indicator Bag, Mask and Star Lumen Oxygen Tubing provides a convenient and cost effective solution.

Through Recovery to Ward

When the laryngeal mask is removed, simply replace the T-Piece with the Aerosol Mask included in the kit for continued oxygen therapy during transfer back to the ward.

		QTY
032-10-640	FCOT T-Piece with Indicator Bag	50
032-10-740	FCOT T-Piece with Indicator Bag, Aerosol Mask and 2.1m Oxygen Tubing	50
032-10-740A	FCOT KIT T-Piece with Indicator Bag, Aerosol Mask and 3m Oxygen Tubing	50
032-10-741	FCOT T-Piece with Indicator Bag, Paediatric Aerosol Mask and 2.1m Oxygen Tubing	50



Post-Operative Oxygen Therapy T-Piece

Fits to the laryngeal mask to deliver 40% oxygen to patients recovering from anaesthesia.

		QTY
032-10-540	40% Oxygen T-Piece with Corrugated Tube	50
032-10-541	Oxygen T-Piece with 24% 28% 31% 35% 40% & 60% Venturi Barrels	50